

Panther Endo Stapler Reload Units

PRODUCT	Endo Articulating reload unit with P3H and P2H technology
PRODUCT CODES	CADD-30EVS, CADD-30ENTS, CADD-45EVS, CADD-45ENTS, CADD-45EPTS, CADD-60EVS, CADD-60ENTS, CADD-60EPTS, CADF-30T, CADF-30N, CADF-45T, CADF-45N, CADF-45S, CADF-45D, CADF-45R, CADF-60T, CADF-60N, CADF-60S, CADF-60D, CADF-60R, CADE-30EVS, CADE-45EVS, CADE-45ENTS
ARTG No	386361

PRODUCT DESCRIPTION

CADD series are articulating reloading units with variable-sized staples, CADE series are articulating reloading units with variable-sized staples and Curved tip. CADF series are articulating reloading units with a fix anvil. The aforementioned reload units can be used with the Panther SMART Endo Stapler and Panther Manual Endo Stapler and delivers Haemostasis and Safety. The enhanced anvil delivers better compression and manipulation of the tissue.

FEATURES:


P3H technology

- Lateral diffusion of fluids
- Stepped cartridge surface
- Better tissue compression
- New knife for each reload unit

P2G technology

- Double gripping technology
- Lateral diffusion of fluids
- Better tissue compression
- Anvil and cartridge with progressive clamping technology for better clamping and less squeezing of the tissue
- New knife for each reload unit

PRODUCT IMAGE



PACKAGING SPECIFICATIONS

• Individually Packed	• Box of 6 UOM	• Sterile Single Use
-----------------------	----------------	----------------------