

Panther SMART Endo Stapler	
PRODUCT	PANTHER SMART Endoscopic motorized cutting stapler, Single-use with Non-Stop articulation system
PRODUCT CODES	PEAS, PEAM, PEAL
ARTG No	386366

PRODUCT DESCRIPTION

The SMART endoscopic motorized cutting stapler, single-use is intended for transection, resection, and/or creation of anastomoses. The instruments have application in multiple open or minimally invasive general, colorectal, bariatric and thoracic surgical procedures.

FEATURES:

- SMART reload identification
- SMART tissue recognition
- SMART firing control
- 15 Firings
- Electronically assisted for precise firing
- Innovative 45 degrees non-stop articulation system
- Universal handle ready to use with most extended portfolio reloads
- Ergonomic universal handle design with 360-degree rotation
- Can be used with P3H and P2G Reload technology
- Sterile, single-use

PRODUCT IMAGE



PACKAGING SPECIFICATIONS

• Individually Packed	• Box of 1 UOM	• Sterile Single Use
-----------------------	----------------	----------------------