

Laparoscopic Monopolar Scissors

PRODUCT	Laparoscopic Monopolar Scissors
PRODUCT CODES	106Y.201 – 33cm, 106Y.202 – 45cm, 106Y.203 -33cm Mini Metz.
ARTG No	365070

PRODUCT DESCRIPTION

The Laparoscopic Monopolar Scissor is designed to cut, dissect and coagulate selected tissue in laparoscopic surgery. It is to be connected through a suitable cable with the monopolar output of an electro-surgical unit/generator.

FEATURES:

- Metzenbaum scissor tips
- Quality, functionality and performance
- Ergonomic handle offers comfort during long periods of use
- Easy to move rotator
- Fully insulated shaft/cannula
- 33cm and 45cm lengths
- Lasting scissor sharpness
- Sterile Single Use
- Competitively priced
- Box of 5

PRODUCT IMAGE



PACKAGING SPECIFICATIONS

• Individually Packed	• Box of 5 UOM	• Sterile Single Use
-----------------------	----------------	----------------------