

Cholangiogram Forceps

PRODUCT	Cholangiogram Forceps – 2 Options
PRODUCT CODE	HD1010116, (11-249G) HD1010116
ARTG No	364393

PRODUCT DESCRIPTION

The cholangiogram forceps are used to grasp the bile ducts, as well as for the placement of a catheter for intra-operative imaging of the gallbladder or the bile ducts. The **Single Use** Cholangiogram Forceps offer guaranteed sterility for every patient, every time.

Features:

- High quality device
- 5mm width x 330mm length
- Ergonomic, comfortable, and light weight design
- Accommodates all catheter sizes up to 5.5FR
- Fined tooth ratchet handle
- Sterile Single Use
- Competitively priced
- Box of 5



PRODUCT IMAGES

HD1010116



(11-249G) HD1010116



PACKAGING SPECIFICATIONS

• Individually Packed	• Box/5 UOM	• Sterile Single Use
-----------------------	-------------	----------------------